

# Dry Granulation Pharma Machine



## Features

- Both the stirring blade and cutter feature variable frequency speed control allowing you to adjust the speed to suit different granulation requirements.
- The PLC-controlled human-machine interface ensures that operations are synchronized with the program, simplifying the overall process through automation.
- The entire process of dry and wet mixing granulation occurs within the same closed container, reducing adhesive usage by 25% and shortening drying time compared to traditional methods.

## Technical Specification

Feeding Motor	0.4 kw
Feeding Speed	0-30 rpm
Whole Grain Motor	0.75kw
Hydraulic Motor	0.75kw
Max oil Pressure	10 Mpa
Weight	450kg
Host Size(L x W x H)	900x760x1650mm

