Rapid Mixing Pharma Granulator



Features

- Equipped with a dripping device, allowing adhesives to be added either by dripping or mixing. This feature is designed to meet the needs of various medicines, saving approximately 20% on adhesive usage while ensuring high-quality granulation.
- The granulator completes the mixing and granulation process in a single step improving efficiency by 4-5 times. The resulting particles are solid and have a high granulation rate.
- The stirring paddle is sealed with a rubber pad, employing a double composite technology that combines oil and air screws to overcome the common issues associated with shaft seal failures.

Technical Specification

>2 Year	
75-210kg	
150-420L	
1150mm	
20-200 rpm	
30kw	
50-3000 rpm	
11kw	
3000KG	
3300mm*2300mm*2550mm	





